

## DeVilbiss IntelliPAP® Specifications

Provided for your reference by  
**cpapXchange.com**

Specification	Standard	Standard Plus	AutoAdjust®	Bilevel	AutoBilevel
<b>Product Number</b>	DV51D	DV53D	DV54D	DV55D	DV57D
<b>Product Dimensions</b>	6.4" H x 6.5" W x 8.4" D	6.4" H x 6.5" W x 8.4" D	6.4" H x 6.5" W x 8.4" D	6.4" H x 6.5" W x 8.4" D	6.4" H x 6.5" W x 8.4" D
<b>Product Weight</b>	4.45 lb with humidifier and PAP	4.45 lb with humidifier and PAP	4.45 lb with humidifier and PAP	4.45 lb with humidifier and PAP	4.45 lb with humidifier and PAP
<b>Pressure Range</b>	3-20 cm H <sub>2</sub> O (.5 increments)	3-20 cm H <sub>2</sub> O (.5 increments)	3-20 cm H <sub>2</sub> O (.5 increments)	3-25 cm H <sub>2</sub> O (.5 increments)	3-25 cm H <sub>2</sub> O (.5 increments)
<b>Universal Power Supply</b>	Yes	Yes	Yes	Yes	Yes
<b>AC Cord Length</b>	8 feet	8 feet	8 feet	8 feet	8 feet
<b>Automatic Altitude Adjustment</b>	Yes	Yes	Yes	Yes	Yes
<b>DC Operable</b>	Yes	Yes	Yes	Yes	Yes
<b>ABOB Analog Output For PSG</b>	Yes (pressure only)	Yes	Yes	Yes	Yes
<b>Automatic Leak Compensation</b>	Yes	Yes	Yes	Yes	Yes
<b>Auto On/Off</b>	Yes (can be enabled or disabled)	Yes (can be enabled or disabled)	Yes (can be enabled or disabled)	Yes (can be enabled or disabled)	Yes (can be enabled or disabled)
<b>Motor Hour Meter</b>	Yes	Yes	Yes	Yes	Yes
<b>Therapy Compliance Meter</b>	Yes	Yes	Yes	Yes	Yes
<b>Standard Warranty</b>	3-year warranty on IntelliPAP	3-year warranty on IntelliPAP	3-year warranty on IntelliPAP	3-year warranty on IntelliPAP	3-year warranty on IntelliPAP
<b>DeVilbiss SmartLink® Compatible</b>	Yes	Yes	Yes	Yes	Yes
<b>SmartCode® Feature</b>	Yes	Yes	Yes	Yes	Yes
<b>Tubing Lengths</b>	6 feet or 10 feet	6 feet or 10 feet	6 feet or 10 feet	6 feet or 10 feet	6 feet or 10 feet



DeVilbiss Healthcare | 100 DeVilbiss Drive | Somerset, PA 15501 | USA  
 800.338.1988 | 814.443.4881 | [www.DeVilbissHealthcare.com](http://www.DeVilbissHealthcare.com)