

Provided for your reference by  
**cpapXchange.com**



## DeVilbiss IntelliPAP® Standard CPAP System

You talked. We listened.  
**Everyone gets what they need.**

The new **DeVilbiss IntelliPAP Standard CPAP System** has been designed with patients and providers in mind, incorporating many of the recommendations solicited through research to **optimize patient adherence**. IntelliPAP's combination of **comfort, education and adherence tracking with SmartCode™** helps to ensure patient compliance – all in a **highly efficient platform** that makes great business sense to you.

From its **reliable electronics and water resistant housing** to its durable Lexan® humidifier chamber and stainless steel heat transfer plate, IntelliPAP is **built for years of dependable performance and patient satisfaction**. Just what you would expect from DeVilbiss – the widely recognized leader in respiratory technology.

DeVilbiss IntelliPAP  
Standard CPAP Unit  
Model DV51D-HH

### Performance Features

- Quietest **CPAP Unit on the market** – ONLY 26 dBA!\*
- Integrated Heated Humidification
- **Automatic Altitude Adjustment** up to 9,000 feet
- Exclusive humidifier chamber design **maximizes humidification and minimizes leaks**
- Exceptionally small footprint, **smallest integrated unit** among top-selling CPAP
- Delay pressure ramp
- Best-in-industry warranty

## Design Features

- Clean, sophisticated design
- Whisper-quiet operation – 26 dBA\*
- Exceptionally small footprint fits easily into the home environment and carry-on luggage
- Discreet carry bag included that can be easily expanded to hold additional accessories
- Integrated heated humidification system easily docks into the bottom of the unit without requiring a separate hose or power cord
- Humidifier reservoir is easy to remove and replace and includes built-in, overflow protection
- Water ingress protection prevents humidifier water from damaging the CPAP – great for travelers
- Patient air tubing extends from back of unit to reduce clutter and enhance visual appeal
- Accommodates various tubing lengths (6 and 10 feet)
- Large, slip-resistant pads on bottom of unit to minimize movement on night table during the night
- Eight-foot power cord length provides greater positioning flexibility – no bulky power brick on cord
- Easy and intuitive set up – three, large buttons control nearly all functions for simplified patient operation
- Large, backlit LCD and easy-to-follow patient menu
- Auto on/off automatically starts the flow of air after breathing into the mask and stops the flow of air when mask is removed
- Delay pressure ramp eases the machine to the prescribed pressure – allowing the user to fall asleep at a lower pressure (0-45 minutes in 5 minute increments and .5 increments for pressure)
- SmartCode™ compliance tracking encodes usage data that can be communicated to a clinician and then translated
- Reminders alert to the need to replace or check parts
- Extra air inlet and fine particle filter included
- Industry leading, 3-year warranty

## Patient Education and Marketing Support

- User video includes an introduction to CPAP therapy with tips for acclimating to, and optimizing, therapy, as well as how to care for and use the IntelliPAP system
- Patient instruction guide demonstrates how easy it is to use IntelliPAP, while providing a comprehensive overview of sleep apnea, its treatment and tips for acclimating to therapy



DeVilbiss IntelliPAP Standard CPAP System includes attractive travel bag

## Technical Specifications

\* Tested per ISO 17510

<b>Product Dimensions</b>	6.4" H x 6.5" W x 8.4" D	<b>Automatic Leak Compensation</b>	Yes
<b>Product Weight</b>	4.45 lbs w/humidifier and PAP	<b>Auto On/Off</b>	Yes (can be enabled or disabled)
<b>Pressure Range</b>	3-20 cm H2O (.5 increments)	<b>Motor Hour Meter</b>	Yes
<b>Universal Power Supply</b>	Yes	<b>Therapy Compliance Meter</b>	Yes
<b>AC Cord Length</b>	8 feet	<b>Warranty</b>	3-year warranty on IntelliPAP 2-year warranty on humidification system
<b>Automatic Altitude Adjustment</b>	Yes	<b>DeVilbiss SmartLink™ Compatible</b>	Yes
<b>DC Operable</b>	Yes	<b>SmartCode™ Feature</b>	Yes
<b>ABOB Analog Output for PSG</b>	Yes (pressure only)		



DeVilbiss Healthcare | 100 DeVilbiss Drive | Somerset, PA 15501  
800.338.1988 | 814.443.4881 | [www.DeVilbissHealthcare.com](http://www.DeVilbissHealthcare.com)