

## HOW TO IDENTIFY HUMIDIFIER CHAMBER

Humidification is an important part of your sleep therapy treatment to aid in comfort of use. Your Respirationics M Series sleep therapy system includes a removable water chamber that should be replaced as instructed by your Homecare Provider. Since the release of the M Series sleep system, Respirationics has enhanced the replaceable water chamber unit to enable easier use and a better seal. The following pictures will help you identify for your homecare provider the specific water chamber model that is compatible with your humidification platform.

With your current chamber in hand, please order from your Homecare Provider, the part number that corresponds with the pictures in one of the three columns.



**BLUE LATCH**



**BLACK LATCH**



**BLACK LATCH**



**BLUE MIDDLE SEAL**




**BLACK MIDDLE SEAL**



**BLACK MIDDLE SEAL**



**SIDE-BY-SIDE  
OPENINGS**



**SIDE-BY-SIDE  
OPENINGS**



**STAGGERED  
OPENINGS**

This chamber has been replaced with a new model as identified in the center column. Please order part number 1003757.

R2B Chamber Kit  
Please order part number  
1003757

R3B Chamber Kit  
Please order part number  
1003758