

RESPIRONICS®

Provided for your reference by  
**cpapXchange.com**



## COMFORT. FOR US, IT'S PERSONAL.

With its StabilitySelector™ and ball-and-socket headgear attachments, the ComfortSelect™ Nasal Mask delivers the kind of secure, personalized fit that your patients demand. Its textured, dual-layer cushion design and adjustable forehead arm help deliver a comfortable fit and guard against mask leaks. It's one more example of how the Comfort Series™ of Masks helps fit 100% of your patients 100% of the time.



ComfortSelect™

## ComfortSelect™

**Adjustable StabilitySelector™ helps make this Comfort Series™ mask a comfortable, leak-free fit for your patients – and your business.**

The ComfortSelect Nasal Mask is an ideal mask for a wide range of patients. Its adjustable StabilitySelector and soft silicone forehead pad help patients find the right fit to reduce pressure points and mask leaks – without overtightening the headgear. And once you find that perfect fit, the ComfortSelect's ball-and-socket headgear attachments make it easier to keep it – they allow you to remove the mask from the headgear without loosening the headgear straps. Once it's on, even gentle tugging will not result in mask leaks. The fit is that secure. And that makes compliance easier to come by.


We've also made sure that the ComfortSelect is a mask that you can feel secure recommending. It offers many features that have been proven popular in other Respironics masks, like its dual-layer mask cushion and silent exhalation port. ComfortSelect uses fewer mask sizes to fit the same population as comparable masks, and it has minimal parts to lose or break, both of which are particularly comforting when you have a crowded stockroom to manage.

To learn more about the ComfortSelect and the entire Comfort Series of Masks, visit [www.comfortseries.respironics.com](http://www.comfortseries.respironics.com).

## Ordering information

Item	Small	Medium	Small/Wide	HCPCS
ComfortSelect Mask with Headgear	1007930	1007919	1007931	A7034, A7035
ComfortSelect Mask Only	1007933	1007932	1007934	A7034
ComfortSelect Mask with Headgear and Extra Cushion	1035966	1035967	1035968	A7034, A7035, A7032
Medium ComfortSelect FitPack™ with Headgear and Small and Small/Wide Cushions		1026880		A7034, A7035

Parts and Accessories	Part #	HCPCS
Small Cushion with Retaining Ring	1007936	A7032
Medium Cushion with Retaining Ring	1007935	A7032
Small/Wide Cushion with Retaining Ring	1007937	A7032
Deluxe Headgear	1002800	A7035
Deluxe Headgear, Small	1018755	A7035
Premium Headgear with EZ Peel™ Tabs	1033678	A7035
Forehead Support and Pad	1018401	
Silicone Forehead Pads (4-pack)	1009722	
Ball-and-Socket Headgear Swivel Clips (4-pack)	1009723	
Exhalation Elbow Kit	1004395	A7045
Mask Swivel	7041	
Oxygen Enrichment Attachment	312710	
Sizing Gauge	1012104	



ComfortSelect uses Respironics' one-size-fits-most **Deluxe Headgear**. This unique headgear with Velcro® tabs easily adjusts in eight places for a secure, custom fit. The Deluxe Headgear evenly distributes tension, adds stability and helps to prevent mask leaks. The ComfortSelect and Headgear disassemble easily for cleaning and replacement. The Headgear can be used with most Respironics masks.



## RESPIRONICS®

Customer Service: 1-800-345-6443 or 724-387-4000  
 Respironics Europe, Africa, Middle East: +33-1-47-52-30-00  
 Respironics Asia Pacific: +852-3194-2280  
[www.comfortseries.respironics.com](http://www.comfortseries.respironics.com)

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician.

Respironics is a registered trademark, and ComfortSelect, Comfort Series, EZ Peel, FitPack and StabilitySelector are trademarks of Respironics, Inc. and its affiliates. Velcro® is a trademark of Velcro Industries. Latex Free. Patent Pending.

©2007 Respironics, Inc. and its affiliates. All rights reserved.

Hoech JJ 4/30/07 MCI 4101095 PN 1044393