

Provided for your reference by

**cpapXchange.com**

RESPIRONICS®



## IT'S A FIRST-CLASS UPGRADE FOR PATIENTS WHO NEED A FULL-FACE MASK.

The ComfortGel™ Full mask — our first full-face gel mask with a replaceable cushion — reflects a commitment to innovative, practical product development that's been a Respiroics hallmark for more than 30 years. Our popular blue gel and sure seal technology (SST™) will make it a patient favorite. Homecare providers also will appreciate its all-new design and our new approach to cushion resupply. We're confident you'll think the ComfortGel Full mask is in a class of its own.



**ComfortGel™ Full**

Envisioning tomorrow. Improving today.™

# ComfortGel™ Full

## Proven comfort and an all-new design make this full-face mask our best ever.

Respironics is leading the way to help more people get the healthy, restful sleep they need. We do that by combining the new with the proven — creating products that set new standards and quickly earn your trust. Our latest full-face mask is a perfect example.

The ComfortGel Full mask combines an all-new design with some very familiar — and popular — components. Our comfortable forehead cushion and enhanced StabilitySelector® combine with our popular blue gel replaceable cushion to create a level of comfort and stability no other full-face mask can approach. The unique sure seal technology (SST™) provides a membrane-thin layer over the gel cushion to follow the contours of the face better, allowing wearers to quickly achieve the comfortable seal they need.

The ComfortGel Full mask excels in a sleep lab setting thanks to its removable gel cushion, superior seal and easy-to-clean design. The DuoPack™ replacement kit offers homecare providers two same-size cushions for easier stocking and improved supply replacement opportunities. Patients needing to wear a full-face mask finally get the comfortable, robust-sealing mask they need for restful sleep. Solutions that improve the lives of patients, sleep professionals and homecare providers — one more way that Respironics is envisioning tomorrow and improving today.

### ORDERING INFORMATION


ITEM	SMALL	MEDIUM	LARGE	HCPCS
ComfortGel Full Mask with Headgear	1040135	1040136	1040137	A7030 A7035
ComfortGel Full Mask	1040568	1040569	1040570	A7030
DuoPack: ComfortGel Full Mask and Headgear plus 1 extra ComfortGel Full Cushion and 1 extra SST Flap	1002854	1002855	1002856	A7030 A7035 A7031

### REPLACEMENT CUSHIONS/ FLAPS

SMALL	MEDIUM	LARGE	HCPCS	
ComfortGel Full Cushion	1040798	1040799	1040800	A7031
SST Flap only	1040795	1040796	1040797	
ComfortGel Full Cushion and SST Flap	1044861	1044862	1044863	A7031

### PARTS & ACCESSORIES

PART NUMBER	HCPCS
ComfortGel Full Premium Headgear	1040138 A7035
Silicone Pillow Forehead Pad	1040114
Entrainment Valve Swivel	1049398
Quick Clips, Blue (2 per pack)	1040139
Port Cap (5 per pack)	1010872
Sizing Gauge	1049507



**RESPIRONICS®**

Envisioning tomorrow. Improving today.™

Customer Service: 800-345-6443 or 724-387-4000  
 Respironics Europe, Africa, Middle East: +33-1-47-52-30-00  
 Respironics Asia Pacific: +852-3194-2280  
[www.respironics.com](http://www.respironics.com)

CAUTION: U.S. federal law restricts these devices to sale by or on the order of a physician.

Respironics and StabilitySelector are registered trademarks; ComfortGel, DuoPack, Quick Clip and SST are trademarks; and *Envisioning Tomorrow. Improving Today.* is a service mark of Respironics, Inc. and its affiliates. Patents Pending.

©2008 Respironics, Inc. and its affiliates. All rights reserved.

Hoehc JJ 12/05/07 MCI 4101436 PN 1049615