


Note: For proper mask sizing print out this page at 100%. Do NOT select the 'Shrink oversized pages to paper size' option.
The box below should measure 1 inch square.

A large diagram showing a side profile of a human face with measurement lines for width (W) and height (L) of the nose and mouth area. Below the face are two smaller diagrams showing the width (W) and height (L) of a mask cushion. To the right of these are four mask cushion shapes labeled L, M, S, and P, each with a horizontal double-headed arrow and the text "cut out along lines" inside. At the bottom center is the Respironics logo, followed by the text "RESPIRONICS® ComfortLite™ Simple Cushion Sizing Gauge" and the website "www.comfortseries.respironics.com".

REF 1018624

RESPIRONICS®

ComfortLite™ Simple Cushion Sizing Gauge

www.comfortseries.respironics.com



Note: For proper mask sizing print out this page at 100%. Do NOT select the 'Shrink oversized pages to paper size' option.
The box below should measure 1 inch square.

