

RESPIRONICS®

Provided for your reference by  
**cpapXchange.com**



## DESIGNED TO BE THE PERFECT MASK FOR PATIENTS WHO HATE MASKS.

Some people simply can't be compliant with therapy that involves a traditional nasal mask because they simply can't get a comfortable fit. We made that insight the inspiration behind new ComfortLite™ 2. We give patients the choice of three different cushion types for a more personal, more comfortable fit. And we use the patient's forehead for stability so it barely touches the face. That, along with improved headgear and stronger, more comfortable components, makes the ComfortLite 2 the ideal choice for many patients. It's one more way that Respironics turns ideas and insights into practical solutions that help you fit 100% of your patients 100% of the time.

ComfortLite™ 2



# ComfortLite™ 2

**More cushion options and no pressure points – designed to help increase patient comfort and compliance.**

The unique design of the new ComfortLite 2 makes it uniquely qualified to satisfy the comfort and compliance needs of patients. Unlike other masks that use the nostrils for seal and stability, the ComfortLite 2 uses its new “Halo” headgear and forehead mechanism for stability. This allows the cushion to provide a comfortable, reliable seal with no pressure points – and that’s why we put special focus on the cushion.

The ComfortLite 2 reaches an all new level of personalized comfort with three cushion styles available: Pillows Cushion, Simple Cushion™ and Direct Seal™ Cushion. Our mask combo packs are available in a variety of styles and sizes.

For improved comfort and durability, we’ve replaced the plastic side straps with a flexible, sleep-friendly fabric. The crown swivel is also made of a new, easy-rotating material, and it’s been moved forward to reduce contact with the patient’s bed pillow. Improvements like this help make the ComfortLite 2 our best lightweight mask ever.

## Ordering Info

Item	Part Number Incl. Headgear	HCPCS	Part Number Mask Only	HCPCS
<b>Mask Subassembly – no cushions</b>	1025225	A7035*	1025194	NA*
<b>COMBO PACKS</b>				
<b>ComfortLite 2 Mask Combo</b> Sizes: Pillows: S, M; Simple: S, M	1030494	A7034 A7035	1030496	A7034
<b>ComfortLite 2 Mask Combo</b> Sizes: Pillows: M, L; Simple: M, L	1030495	A7034 A7035	1030497	A7034
<b>ComfortLite 2 Mask Combo</b> Sizes: Direct Seal: 4, 5; Simple: S, M	1025112	A7034 A7035	1025121	A7034
<b>ComfortLite 2 Mask Combo</b> Sizes: Direct Seal: 5, 6; Simple: M, L	1025113	A7034 A7035	1025122	A7034
<b>MASK PACKS WITH ONE CUSHION STYLE</b>				
<b>ComfortLite 2 with Pillows Cushion only</b> Sizes: S, M, L	1030493	A7034 A7035	1030498	A7034
<b>ComfortLite 2 with Simple Cushion only</b> Sizes: S, M, L	1025111	A7034 A7035	1025120	A7034
<b>ComfortLite 2 with Direct Seal Cushion only</b> Sizes: 4, 5, 6	1025127	A7034 A7035	1025123	A7034

\*Requires addition of mask cushion for HCPCS-A7034

Replacement Parts	Part Number	HCPCS
Halo Headgear (including strap)	1025117	A7035
Halo Headgear (strap only)	1025114	NA
Pillows Cushion – Petite (1/pack)	1030502	A7032
Pillows Cushion – Small (1/pack)	1030503	A7032
Pillows Cushion – Medium (1/pack)	1030504	A7032
Pillows Cushion – Large (1/pack)	1030505	A7032
Simple Cushion – Petite (1/pack)**	1021835	A7032
Simple Cushion – Small (1/pack)**	1021830	A7032
Simple Cushion – Medium (1/pack)**	1021823	A7032
Simple Cushion – Large (1/pack)**	1021821	A7032
Direct Seal Cushion – Size 1 (1/pack)	1021802	A7032
Direct Seal Cushion – Size 2 (1/pack)	1021806	A7032
Direct Seal Cushion – Size 3 (1/pack)**	1021807	A7032
Direct Seal Cushion – Size 4 (1/pack)**	1021815	A7032
Direct Seal Cushion – Size 5 (1/pack)**	1021816	A7032
Direct Seal Cushion – Size 6 (1/pack)**	1021817	A7032
Pillows Cushion Sizing Gauge (5/pack)	1032007	NA
Pillows Cushion Sizing Gauge (20/pack)	1031154	NA
Simple Cushion Elbow (5/pack)	1014206	NA
Simple Cushion Sizing Gauge	1018624	NA
Direct Seal Cushion Sizing Gauge (20/pack)	1018665	NA
Small Accessory Pouch	1027999	NA
ComfortLite 2 Fitting Guide	1029161	NA

\*\* Also available in 5/packs.

Updated 360° crown swivel is made of a new material to swivel more easily and move with the patient.

Tubing attachment point moved forward toward forehead, minimizing contact with pillow – and maximizing stability.

Multiple position angle adjuster and rotating vertical adjuster provide optimal cushion placement for personalized seal and superior comfort.

Headgear straps are now made with a softer, sleep-friendly fabric to eliminate pressure points and improve comfort for all sleepers.

Flexible adjustment band on Pillows Cushion and Direct Seal™ optimizes cushion angle for the best seal.



A more traditional cushion option, the Simple Cushion™ has a freely rotating ball-and-socket elbow/cushion assembly, which enables proper cushion tilt for a personalized seal.

## Three Nasal Interface Options



NEW Pillows Cushion inserts into the nostrils for a quick seal and easy patient set-up.



Simple Cushion™ is a tip-of-the-nose nasal cushion that is smaller to reduce bulk and pressure points.



Direct Seal™ Cushion offers a comfortable seal without insertion into the nostrils.

# RESPIRONICS®

Customer Service: 800-345-6443 or 724-387-4000

Respironics Europe: +33-1-47-52-30-00

Respironics Asia Pacific: +81-3-5280-9611

www.comfortseries.respironics.com

CAUTION: U.S. federal law restricts these devices to sale by or on the order of a physician.

Comfort Series, ComfortLite, Direct Seal, Simple Cushion and Respironics are trademarks of Respironics, Inc. and its affiliates. Patent Pending.

©2006 Respironics, Inc. and its affiliates. All rights reserved.

Hoech – BK 02/24/06 MCI 4100333

PN 1028746