

## DeVilbiss FlexSet® Nasal CPAP Mask

### One Mask, So Many Options

With so many nasal CPAP mask options on the market, it is cost prohibitive to stock all of them. The DeVilbiss FlexSet mask provides four different cushion options, thus you can minimize the array of masks you stock and **build the mask to suit your patient**. Additionally, the tubing can be worn straight down or over-the-head for added flexibility.

The FlexSet is amazingly intuitive to use and care for. This innovative interface features a unique forehead adjustment and pivoting design that help to ensure a **secure, leak-free fit on most facial profiles**. Furthermore, its lightweight design and minimal headgear make the mask **extremely comfortable** for the patient to wear.

At DeVilbiss Healthcare it is our goal to provide you with options that will satisfy your needs as well as the needs of your patients.



DeVilbiss FlexSet  
Nasal CPAP Mask  
9354D

### Performance Features

- Patented **adjustable exhaust port** allows patient to direct escaping air away from their chest or bed partner
- One-piece, **self-anchoring headgear** is easy to use
- Dual-pivoting forehead rest and mask bridge for greater comfort and improved fit
- Two ultra-soft Breath-O-Prene® forehead pads included for optimal patient comfort
- Hose connection swivel elbow that **rotates 360 degrees** or can be locked in an upright position to run tubing over the head.

## Design Features

- 100% latex free
- Soft, replaceable cushion provides a comfortable low pressure seal. Available in standard and shallow sizes and silicone or gel material.
- Quick disconnect for ease of use – tubing can be disconnected from the mask swivel.
- Minimal headgear offers a secure fit and features a quick release buckle.
- Supplemental oxygen port.
- The FlexSet mask provides a modular and easy-to-use nasal CPAP interface that allows users to build a mask to suit their needs.
- There are currently 4 different cushion options available (see accessory and replacement parts listing below) and each FlexSet mask comes complete with 2 different forehead pad options.
- The FlexSet Mask is dishwasher safe (with cushion, headgear and forehead pad removed).

## Reimbursement information

Product Description	Code Description	HCPCS Code	Frequency	Allowable Range
FlexSet Mask	Nasal Application Device (mask or cannula type), used with positive airway pressure device, with or without head strap	A7034	1 per 3 months	\$99.00 - \$117.64 initial set-up
Replacement Cushion	Replacement cushion for Nasal Application Device, each	A7032	2 per 1 month	\$34.45 - \$40.53
Replacement Headgear	Headgear	A7035	1 per 6 months	\$33.79 - \$39.75
Chin Strap	Chinstrap	A7036	1 per 6 months	\$15.47 - \$18.20

## Accessory and Replacement Parts

ITEM	DEVILBISS PART #
FlexSet Mask Standard	9354D
FlexSet Mask Shallow	9354S
FlexSet Mask Gel	9354G
FlexSet Mask Gel Shallow	9354GS
FlexSet Silicone Standard Cushion	9352D-601
FlexSet Gel Standard Cushion	9352G-601
FlexSet Silicone Shallow Cushion	9352S-601
FlexSet Gel Shallow Cushion	9352GS-601
FlexSet Headgear	9353D-619
FlexSet Chin Strap	7351D-647
FlexSet Pad/Vent	9353D-604



Standard and Shallow Gel Cushions  
9352G-601, 9352GS-601



DeVilbiss Healthcare | 100 DeVilbiss Drive | Somerset, PA 15501 USA  
800.338.1988 | 814.443.4881 | [www.DeVilbissHealthcare.com](http://www.DeVilbissHealthcare.com)