

Provided for your reference by  
**cpapXchange.com**



# *Transforming sleep.*

## **Forma™**

### **FULL FACE MASK**

The Forma™ Full Face Mask forms a superior seal and provides even more comfort with a greater range of movement through the active contouring of its new FlexiFoam™ Cushion. The T-Piece adds stability to the stylized Mask Base, and works in harmony with the highly contoured Silicone Seal and Under-Chin Design to enhance sleep performance. With its ease of use and superior conforming seal, the Forma™ sets a whole new standard.


Transform your patient's sleep experience today, with the Fisher & Paykel Healthcare Forma™ Full Face Mask.

**FlexiFoam™ Cushion & T-Piece Mask Base**




**Fisher & Paykel**  
HEALTHCARE

# Full Face Delivery – *Forma™*



The *Forma™* Full Face Mask offers an improved seal and a more comfortable fit due to the interaction between the new *FlexiFoam™* Cushion and the Under-Chin Design. This alternative to nasal or oral delivery comes fully assembled and ready for use.

## UNIQUE MASK FEATURES FOR GREATER FREEDOM, COMFORT AND COMPLIANCE




### FlexiFit™ Technology Auto-contouring fit

The *FlexiFoam™* Cushion using *FlexiFit™* Technology automatically adjusts to a wide range of nasal and facial contours with simple one-step ease of fitting and no complicated adjustments.




### FlexiFoam™ Cushion Active-contouring cushion

The *FlexiFoam™* Cushion provides greater comfort with the unique design of the soft, light and supportive foam. The *FlexiFoam™* features the advanced active contouring that conforms naturally to your patient's face.




### Glider™ Strap with Easy-Release Cord

The sliding *Glider™* Strap with *Easy-Release Cord* provides freedom of movement and one-step removal. The mask force remains equal when the head rotates from side to side, minimizing the occurrence of leaks. Simply pull the *Easy-Release Cord* for one-step removal.




### Under-Chin Design with Crown Strap

An Under-Chin fit provides extra stability and maintains the seal, even with movement of the jaw during the night. The *Crown Strap* works in conjunction with the Under-Chin Design to eliminate the need for a chin strap.




### T-Piece Mask Base Optimal seal

With no complicated adjustments, the *T-Piece* provides added stability and complements the Under-Chin Design to allow *FlexiFit™* to auto-adjust for hassle-free fitting.



### Advanced Air Diffuser Quiet usage

**Maintenance Free.** This efficient *Air Diffuser* system prevents air from blowing onto the sleeping partner and decreases noise that could potentially disrupt sleep.



### Frosted Silicone Seal Optimizing seal

The enhanced seal technology has a soft, smooth feel and slides easily onto the face. Its frosted texture optimizes mask seal performance.

## CONTENTS

- Complete mask with headgear assembled and ready for use
- Four sizes: Small, Medium and Large (two sizes included) and X-Large, sold separately
- User-friendly Fitting Guide

## REORDER CODE:

400470 Small  
400471 Medium and Large  
400473 X-Large

FOR MORE INFORMATION, PLEASE CONTACT YOUR LOCAL FISHER & PAYKEL HEALTHCARE REPRESENTATIVE:

**Manufacturer** Fisher & Paykel Healthcare Limited, 15 Maurice Paykel Place, East Tamaki, Auckland 2013, PO Box 14 348, Panmure, Auckland 1741, New Zealand  
Tel: +64 9 574 0100 Fax: +64 9 574 0158 Email: info@fphcare.com Web: www.fphcare.com

**Australia** Tel: +61 3 9879 5022 Fax: +61 3 9879 5232 **Austria** Tel: 0800 29 31 23 Fax: 0800 29 31 22 **Benelux** Tel: +31 40 216 3555 Fax: +31 40 216 3554 **China** Tel: +86 20 3205 3486 Fax: +86 20 3205 2132 **France** Tel: +33 1 6446 5201 Fax: +33 1 6446 5221 **Germany** Tel: +49 7182 93777 0 Fax: +49 7182 93777 99 **India** Tel: +91 80 4123 6041 Fax: +91 80 4123 6044 **Irish Republic** Tel: 1800 409 011 **Italy** Tel: +39 06 7839 2939 Fax: +39 06 7814 7709 **Japan** Tel: +81 3 3661 7205 Fax: +81 3 3661 7206 **Northern Ireland** Tel: 0800 132 189 **Spain** Tel: +34 902 013 346 Fax: +34 902 013 379 **Sweden** Tel: +46 8 564 76 680 Fax: +46 8 36 63 10 **Switzerland** Tel: 0800 83 47 63 Fax: 0800 83 47 54 **UK (EU Authorised Representative)** Fisher & Paykel Healthcare UK Ltd, Unit 16, Cordwallis Park, Clivemont Road, Maidenhead, Berkshire SL6 7BU, UK Tel: +44 1628 626 136 Fax: +44 1628 626 146 **USA** Tel: 1 800 446 3908 or +1 949 453 4000 Fax: +1 949 453 4001