



V2 Full Face
CPAP Mask

V2™ Full-Face CPAP Mask

Redesigned for Improved Comfort and Seal

The V2 full face CPAP mask has been **redesigned to be ultimately more comfortable and to provide a stable seal** for a wider array of facial profiles. These improvements include incorporating a membrane inside the cushion for an **improved seal**. In addition, the cushion surface features a series of flanges around the base and is **constructed of a thinner, softer and more comfortable silicone material**. Two structural braces have been added to aid in cushion stability. Finally, **the CO₂ vents were moved** from the bridge of the nose to just above the anti-asphyxia valve, **pointing the airflow perpendicular to the body**.

Other improvements include **no forehead support, single-molded mask frame and cushion that are made of easy contouring silicone rubber**. The one piece design and soft material allows for greater flexibility through a broad range of 5 different mask sizes.

Performance Features

- Headgear uses **new quick-release strap clips** and a foam material for improved comfort
- **360-degree mask swivel** for expanded patient freedom
- Detachable swivel allows user to **disconnect from equipment without removing mask**
- **5 sizes of mask** are available (mask sizes tend to run large)
- **3 sizes of headgear** are available (small, medium, large)

Design Features

- Lighter weight than original design
- Single-molded mask frame and cushion with **integrated chin support**
- Anatomically contoured design with ribbed support to improve seal and comfort
- Integrated facial-sealing flange for improved fit
- **Anti-Asphyxia Valve** that opens at mask pressure of less than 3 cm H₂O
- Swivel port micro-vent holes that dissipate CO₂
- **Detachable swivel port** that allows the hose to be easily removed
- **Easily adjustable 4-strap headgear** made of durable and comfortable polyurethane foam
- Reusable mask is steam autoclaveable and dishwasher safe. Can also be cleaned via cold chemical and pasteurization processes
- Latex-free
- Multi-patient, multi-use products

V2 Ordering Information

(includes headgear)

| PRODUCT DESCRIPTION | PART NUMBER |
|---------------------|-------------|
| Large | 113485 |
| Medium | 113486 |
| Small | 113487 |
| Extra Small | 113488 |
| Petite | 113489 |

Reimbursement Table

| PRODUCT DESCRIPTION | CODE DESCRIPTION | HCPCS CODING | FREQUENCY | ALLOWABLE RANGE |
|----------------------------|---|--------------|----------------|---------------------|
| Face Piece and Swivel Port | Full face mask used with PAP device, each | A7030 | 1 per 3 months | \$148.44 – \$174.64 |
| Headgear | Headgear used with PAP device, each | A7035 | 1 per 6 months | \$31.27 – \$36.79 |
| Face Piece | Face mask interface, replacement for full face mask, each | A7031 | 1 per 1 month | \$54.90 – \$64.59 |

The noted allowable rates are based off of the January 2012 HCPCS Fee Schedule. All information presented is for reference purposes only.

Accessories and Parts

| PRODUCT DESCRIPTION | PART NUMBER | PRODUCT DESCRIPTION | PART NUMBER |
|---------------------|-------------|------------------------------------|-------------|
| Face Pieces | | Face Piece with Swivel Port | |
| Large | 669180 | Large | 113570 |
| Medium | 669181 | Medium | 113571 |
| Small | 669182 | Small | 113572 |
| Extra Small | 669183 | Extra Small | 113573 |
| Petite | 669184 | Petite | 113574 |
| Headgear | | Mask Sizing Gauge | |
| Large | 201494 | | 691143 |
| Medium | 201493 | | |
| Small | 201492 | | |



DeVilbiss Healthcare LLC | 100 DeVilbiss Drive | Somerset, PA 15501 USA
800.338.1988 | 814.443.4881 | www.DeVilbissHealthcare.com