

## Mirage Activa LT

The Mirage Activa LT nasal mask combines personalized fit and comfort with a new compact ActiveCell cushion to accommodate changing sleeping positions and pressures throughout the night, making it the only mask you'll ever need.



Provided for your reference by  
**cpapXchange.com**


**RESMED**

## Mirage Activa™ LT

NASAL MASK

Fits right for an easy night

Visit [www.myresmed.com](http://www.myresmed.com) for more information about the Mirage Activa LT and other ResMed products, as well as fitting templates and other therapy tips.



1 Testing per ISO 3744:1994 Acoustics determination of sound power levels of noise using pressure at 10 cm H<sub>2</sub>O. Quoted percentage comparisons are calculated by converting sound power values from a logarithmic scale to a linear scale

ResMed Corp Poway, CA, USA +1 858 746 2400 or 1 800 424 0737 (toll free), ResMed Ltd Bella Vista, NSW, Australia +61 (2) 8884 1000 or 1 800 658 189 (toll free). Offices in Austria, Belgium, Brazil, China, Finland, France, Germany, Greece, Hong Kong, India, Ireland, Italy, Japan, Malaysia, Mexico, Netherlands, New Zealand, Norway, Portugal, Singapore, South Africa, Spain, Sweden, Switzerland, United Kingdom (see website for details). Activa, ActiveCell and Mirage are trademarks of ResMed Ltd. and are registered in U.S. Patent and Trademark Office. ©2008 ResMed. 1013304/1 08 11



## Personalized fit and comfort

- **New compact ActiveCell™ technology**  
Mirage™ dual-wall cushion with ActiveCell expands and contracts to accommodate more sleeping positions and varying pressures
- **Gentle seal**  
Lightest sealing pressure for a gentle fit that reduces facial marks and eye irritation, and promotes long-term comfort
- **Self adjusting**  
Decouples movement from mask to the face for an uninterrupted night's sleep
- **First-time fit**  
MicroFit dial fine-tunes fit to suit your individual profile for a comfortable seal



ActiveCell expands ...



... and contracts



## Quiet and streamlined

- **Quiet performance**  
49% quieter than Mirage Activa<sup>1</sup>
- **Softly diffused**  
Unique design diffuses air quietly and gently for improved comfort for you and your bed partner
- **Clear field of vision**  
Slender frame design and compact cushion are unobtrusive and feel light on the face

## Intuitively easy

- **Hassle-free seal**  
ActiveCell automatically seals and adjusts, reducing fitting time
- **One-handed adjustment**  
You can easily adjust the MicroFit dial with one hand, even while lying down
- **Easy maintenance of mask**  
Fewer mask parts make it easy to clean, assemble and disassemble
- **Comprehensive user support package**  
Includes travel bag and user guide with mask fitting and care instructions



The MicroFit dial provides 24 positions of comfort at your fingertips