

Technology Highlight: **Cushion and Elbow**

"Very comfortable;
cushion feels like suede."

"I like the fact that the
air [from vents] isn't direct."

Reference: ResMed user trials



The dual-wall, Spring Air cushion is coated with soft-touch frosting to provide comfort and even pressure distribution. Incorporating a softer texture at the nasal bridge—and the absence of an undercushion—reduces soreness, irritation and discomfort commonly felt in this sensitive area. The absence of a cushion clip reduces the total number of parts and simplifies the assembly/disassembly process.



Squeeze-tab elbow

makes it easy to attach and detach the elbow from the frame with just one hand.

Tubing can be attached to the elbow prior to the elbow being attached to the frame. This eliminates the hard snap of a swivel and pushing into the face.

The innovative elbow to cushion assembly is another way Mirage™ FX achieves both sophisticated design and easy functionality.



In-moulded vent design

moulded into the cushion to disperse air gently and quietly away from bed partner.

Sound Power: 26 dbA
(A-weighted sound power level of the mask, internal ResMed testing)

All this in One cushion size that caters for a wide range of facial characteristics.



Effortless and
effective

Mirage™ FX
NASAL MASK

www.resmedmiragefx.com