

Mirage Quattro™ FULL FACE MASK


FOURTH GENERATION TECHNOLOGY

Wide fit range: fourth generation full face mask from ResMed fits more than 95% of users.*

Mirage™ dual-wall cushion with spring air technology: unique cushion accommodates movement during sleep and provides an effective solution for jaw drop.

MicroFit dial: new forehead support adjustment offers personalized fitting and comfort.

Light, quiet and easy to use: these benefits are optimized with streamlined mask design, improved vent design and a range of user-friendly features.





TECHNOLOGY HIGHLIGHT

MicroFit dial

A choice of 24 positions at your fingertips—providing a comfortable, effective fit for over 95% of users*

- **Adjustable:** can be adjusted while on the user's face
- **Precise:** each turn of the dial adjusts the forehead support in small increments
- **Progressive:** moves in towards the face, allowing the user to sense when the fit is most comfortable for them
- **Works in combination with the cushion:** as the dial adjusts the angle of the mask to the face, the cushion adapts to enhance the seal where it's needed most—without the need to overtighten the headgear

USER BENEFITS

Personalized comfort

- MicroFit dial adjusts the mask to suit user's profile
- Dual-wall cushion with spring air technology adapts to facial contours and accommodates movement during sleep (including jaw drop)

Light and quiet

- Streamlined design feels light on the face
- Improved vent design reduces noise

Making it easier

- MicroFit dial and spring air cushion technology help achieve good mask fit—first time, every time
- Fewer parts for easy handling, cleaning, assembly and disassembly

PROFESSIONAL BENEFITS

Fit range increased to over 95% of users*

- MicroFit dial and spring air cushion technology adjust the mask to suit more facial structures
- New XS mask size

First time fit

- First-time fit reduces callbacks and mask refitting
- Improved cushion technology seals quickly
- Simplified parts for quick and easy assembly and disassembly

Reduced inventory

- 33% fewer product codes than Ultra Mirage™ Full Face Mask
- Fewer mask parts

Comprehensive user support package

- Travel and storage bag
- Mask care DVD
- Quick-fitting guide
- User guide



PRODUCT CODES

Mirage Quattro Full Face Mask

X-Small	61200
Small	61201
Medium	61202
Large	61203

REIMBURSEMENT CODES

Mirage Quattro Full Face Mask

HCPSC Code	Descriptor
A7030	Full face mask used with positive airway pressure (PAP) device, 1 per 3 months
A7031	Replacement cushion for full face mask, 1 per month
A7035	Headgear used with positive airway pressure (PAP) device, 1 per 6 months

* Source: internal ResMed anthropometric data