


 Sleepnet®  
Fine, handcrafted sleep masks

**Mojo™**

**Full Face Mask  
Size Gauge**

Small 50835  
Medium 50836  
Large 50837



Example: Patient takes a SMALL

Hold size gauge up to patient's face and measure from indent on nose bridge to indent under lips. Relax jaw, with mouth slack and slightly open (about 1 1/2.5cm).

Provided for your reference by  
**cpapXchange.com**

90102 Rev. B



Print out template at 100%.  
Make sure not to choose "shrink to fit" when printing.  
To ensure your template is the correct size, measure the "1 Inch" square above.