


Nasal-Aire[®] II Interface

Something New is in the Aire...

The Nasal-Aire II is the **next generation of the ever-popular Nasal-Aire** interface. The Nasal-Aire II shares many of the **comfort and compliance benefits** of the original Nasal-Aire, adding improved performance and a reduction in sound level.

This lightweight interface is comfortable for patients, simple to use and easy to care for. Nasal-Aire II is **sold with all six sizes for easy fitting** and includes headgear for stability.



Nasal-Aire II
Nasal CPAP Mask
K2A

Performance Features

- Ergonomically designed to fit any face with no leaks
- Virtually eliminates uncomfortable pressure points on the face allowing the interface to be worn for any duration without discomfort
- Secure fit design allows the user to sleep in any position
- Minimal design reduces feelings of claustrophobia
- New ultra-soft nasal prong material for improved patient comfort

Design Features

- 100% latex free
- Medium or large size fits ninety percent of patients
- Can be worn over the ears or down the chest
- Small and quiet with virtually no noise
- No nasal bridge or forehead contact eliminates nasal bridge sores and air leaks near eyes
- Swivel tubing coupler to accommodate free movement and optimal performance during the night
- Supplemental oxygen ports

Reimbursement Information

Product Description	Code Description	HCPCS Coding	Frequency	Allowable Range
NASAL-AIRE II INSERT	Cushion for nasal application device	A7032	Two per month	\$34 -\$40
NASAL-AIRE II APPLICATION	Nasal interface (mask or cannula type) used with positive airway pressure device, with or without head strap	A7034	One per 3 months	\$99 -\$117
HEADGEAR	Headgear used with positive airway pressure device	A7035	One per 6 months	\$33 -\$39

Accessories and Replacement Parts

PART NUMBER	PRODUCT DESCRIPTION	PART NUMBER	PRODUCT DESCRIPTION
K2A	Nasal-Aire II All Size Kit, Includes XS, SM, MD, MP LG, XL*	XS01	Nasal-Aire II Cannula Only, XS
K2LG	Nasal-Aire II New Patient Setup Kit, MP, LG, XL*	SM02	Nasal-Aire II Cannula Only, SM
K2MD	Nasal-Aire II New Patient Setup Kit, SM, MD, MP*	MD03	Nasal-Aire II Cannula Only, MD
K2MP	Nasal-Aire II New Patient Setup Kit, MD, MP, LG*	MP04	Nasal-Aire II Cannula Only, MP
K2SM	Nasal-Aire II New Patient Setup Kit, XS, SM, MD*	LG05	Nasal-Aire II Cannula Only, LG
PA401	Nasal-Aire II Petite, Size A	XL06	Nasal-Aire II Cannula Only, XL
PB402	Nasal-Aire II Petite, Size B	NAT320	Nasal-Aire II Frame Kit
PC403	Nasal-Aire II Petite, Size C	X2105	Nasal-Aire II Headgear with OHS Assembly
PD404	Nasal-Aire II Petite, Size D	P430	Nasal-Aire II Petite – Headgear
PE405A	Nasal-Aire II Petite, Size E	PA411	Nasal-Aire II Petite – Cannula Only, Size A
XS301	Nasal-Aire II, XS**	PB412	Nasal-Aire II Petite – Cannula Only, Size B
SM302	Nasal-Aire II, SM**	PC413	Nasal-Aire II Petite – Cannula Only, Size C
MD303	Nasal-Aire II, MD**	PD414	Nasal-Aire II Petite – Cannula Only, Size D
MP304	Nasal-Aire II, MP**	PE415A	Nasal-Aire II Petite – Cannula Only, Size E
LG305	Nasal-Aire II, LG**	X2007	Sizing Gauge – Nasal-Aire II
XL306A	Nasal-Aire II, XL**	X4007	Sizing Gauge – Nasal-Aire II Petite

* Includes headgear, frame, bolo and multiple nasal cannulas.

** Includes headgear, frame, nasal cannula and bolo.



DeVilbiss Healthcare | 100 DeVilbiss Drive | Somerset, PA 15501 USA
800.338.1988 | 814.443.4881 | www.DeVilbissHealthcare.com