



### Determining your best fit for the Quattro™ FX / Quattro™ FX For Her Mask

Measure from a point at the top of the nasal bridge down to the indentation below your lips and above your chin as shown in the illustration above.

Use the measurements below as a guide to find your best size.

Please note that every user is different. These sizes should be used as a guide only.

<b>Large Quattro FX</b>	<b>4 3/8 inches</b>	<b>111 millimeters</b>
<b>Medium Quattro FX / Quattro FX For Her</b>	<b>4 inches</b>	<b>101 millimeters</b>
<b>Small Quattro FX / Quattro FX For Her</b>	<b>3 9/16 inches</b>	<b>90 millimeters</b>