

Mask and accessories materials list

Includes all Lab versions except where noted.


All parts are **latex free**

Product	Material
ComfortGel Blue	
Faceplate	Polycarbonate
Gel cushion	Urethane gel/EVA
Gel cushion outer shell	Urethane film
Outer flap	Silicone
Forehead cushion	Silicone
Exhalation elbow and swivel	Polycarbonate
Elbow split washer	Polycarbonate
Forehead support arm and bracket	Polycarbonate
Headgear clips	Hytrel
Headgear (premium)	UBL/Urethane foam/Lycra

Product	Material
ComfortGel	
Faceplate	Polycarbonate
Gel cushion	Urethane gel
Gel cushion outer shell	Urethane film
Outer flap	Silicone
Forehead cushion	Silicone
Exhalation elbow and swivel	Polycarbonate
Elbow split washer	Polycarbonate
Forehead support arm and bracket	Polycarbonate
Headgear clips	Hytrel
Headgear (premium)	UBL/Urethane foam/Lycra

Product	Material
EasyLife	
Faceplate	Polycarbonate
Inner seal cushion	Silicone
Outer support cushion	Silicone
Diaphragm cushion	Silicone
Exhalation elbow	Polycarbonate
Standard elbow (SE)	Polycarbonate
Elbow split washer	Polycarbonate
Headgear	UBL/Urethane foam/Lycra

Product	Material
ComfortGel Full	
Faceplate	Polycarbonate
Gel cushion	Urethane gel
Gel cushion outer shell	Urethane
Outer flap (SST)	Silicone
Cushion flap ring	Polycarbonate
Forehead cushion	Silicone
Exhalation elbow and swivel	Polycarbonate
Elbow split washer	Silicone
Entrainment valve flapper	Silicone
Forehead support arm and bracket	Polycarbonate
Headgear clips	Hytrel
Headgear (premium)	UBL/Urethane foam/Lycra



Product	Material
FullLife	
Frame	Hytrel
Outer flap (SST)	Silicone
Entrainment elbow and swivel	Polycarbonate
Elbow split washer	Polycarbonate
Entrainment valve flapper	Silicone
Headgear	UBL/Urethane foam/Lycra and Nylon

Product	Material
FitLife	
Faceplate	Polycarbonate
Cushion	Silicone
Exhalation elbow and swivel	Polycarbonate
Faceplate valve hub	Polycarbonate
Accessory swivel	Polycarbonate
Elbow split washer	Polycarbonate
Entrainment valve flapper	Silicone
Port cap	Polycarbonate
Headgear clips	Acetal
Headgear	UBL/Urethane foam/Lycra
Headgear (lab version)	Laminate composition of nylon and polyester polyurethane foam

Product	Material
OptiLife	
Frame	Hytrel
Cushion (pillows and CradleCushion)	Silicone
Tubing	Ethylene vinyl acetate (EVA)
Swivel	Polycarbonate
Headgear	UBL/Urethane foam/Lycra

Product	Material
ComfortFusion	
Faceplate	Polycarbonate
Cushion	Silicone
Forehead cushion	Silicone
Exhalation elbow and swivel	Polycarbonate
Elbow split washer	Polycarbonate and/or PETG
Forehead support arm and bracket	Polycarbonate
Headgear clips	Hytrel
Headgear	UBL/Urethane foam/Lycra

Product	Material
ComfortLite 2	
Cushion (pillows, Simple Cushion, Direct Seal)	Silicone
Simple Cushion elbow	Polycarbonate
Direct Seal cushion cradle	Polycarbonate
Direct Seal cushion support band	Aluminum
Crown swivel	Polycarbonate
Flex tube	EVA copolymer
Pivot tube insert	Polycarbonate
Vertical adjustment nut	Polypropylene
Horizontal adjustment tube	Polypropylene
Pivot mount strap	Polycarbonate
Pivot tube	Polypropylene
Headgear halo	Nylon/polyester
Headgear padding	Nylon/Lycra laminated foam
Headgear strap	Nylon/Lycra laminated foam

Product	Material
ComfortFull 2	
Faceplate	Polycarbonate
Cushion	Silicone
Forehead cushion	Silicone
Exhalation elbow and swivel	Polycarbonate
Entrainment valve flapper	Silicone
Port cap	PVC
Headgear (premium)	UBL/Urethane foam/Lycra
Headgear clips	Hytrel

Product	Material
ComfortSelect	
Faceplate	Polycarbonate
Cushion	Silicone
Forehead cushion	Silicone
Exhalation elbow and swivel	Polycarbonate
Elbow split washer	Polycarbonate
Headgear clips	Hytrel
Headgear (deluxe)	Velstretch/Lycra laminate

ComfortClassic	
Faceplate	Polycarbonate
Cushion	Silicone
Gel forehead pad	Urethane gel encased in a polypropylene sheet
Exhalation elbow and swivel	Polycarbonate
Elbow split washer	Polycarbonate and/or PETG
Port cap	Flexible vinyl (PVC)
Headgear (deluxe)	Velstretch/Lycra laminate

Respironics Profile Lite, Profile Lite Youth, Small Child Profile Lite	
Faceplate	Polycarbonate
Cushion outer shell	Urethane
Port cap	Flexible vinyl (PVC)
Exhalation elbow and swivel	Polycarbonate
Spacer	White polyethyl foam
Headgear (deluxe)	Velstretch/Lycra laminate
Child SoftCap headgear	Velstretch/Lycra laminate
Headgear (RS)	UBL/Urethane foam/Lycra

Respironics Simplicity nasal mask	
Faceplate	Polycarbonate
Cushion	Silicone
Tubing	Ethylene vinyl acetate (EVA)
Deluxe softcap	Velstretch/Lycra nylon
Headgear	Neoprene/nylon

Respironics Total face mask	
Faceplate	Polycarbonate
Cushion	Silicone
Seal ring	Thermoplastic polyester elastomer
Valve body	Polycarbonate
Headgear	Neoprene/nylon

Product	Material
Lab masks (LabClassic, LabFull, LabGel, LabSelect)	
<i>Materials as listed above for reusable masks, except:</i>	
Foam forehead pad	Polyurethane foam
Disposable headgear	Polyurethane foam lined with nylon cloth (40 denier tricut)

Headgear	
Headgear (premium)	UBL/Urethane foam/Lycra
Headgear (RS)	UBL/Urethane foam/Lycra
Headgear (deluxe)	Velstretch/Lycra laminate
SimpleStrap:	
Mesh	Polyester
Strap loop/hook	Velstretch/nylon hook
Softcap:	
Cap	Lycra
Strap loop/hook	Velstretch/nylon hook

Tubing	
Standard gray tubing:	
Cuffs	Silicone rubber or equivalent thermoplastic elastomer
Wall	Polyester elastomer or ethylene octene copolymer
Helix	Polyester elastomer or LDPE or HDPE or equivalent
Performance tubing:	
Cuffs	Thermoplastic elastomer
Wall	Ethylene octene copolymer
Helix	Polypropylene

Philips Healthcare is part of
Royal Philips Electronics

How to reach us

www.philips.com/healthcare
healthcare@philips.com
fax: +31 40 27 64 887

Asia
+852 2821 5888

Europe, Middle East, Africa
+49 7031 463 2254

Latin America
+55 11 2125 0744

North America
+1 425 487 7000
800 285 5585 (toll free, US only)

Philips Respironics
1010 Murry Ridge Lane
Murrysville, PA 15668

Customer Service
+1 724 387 4000
800 345 6443 (toll free, US only)

Respironics Asia Pacific
+65 6882 5282

Respironics Australia
+61 (2) 9666 4444

Respironics Deutschland
+49 8152 93 06 0

Respironics Europe, Middle East, Africa
+33 1 47 52 30 00

Respironics France
+33 2 51 89 36 00

Respironics Italy
+39 039.203.1

Respironics Sweden
+46 8 120 45 900

Respironics Switzerland
+41 6 27 45 17 50

Respironics United Kingdom
+44 800 1300 845

www.philips.com/respironics

Please visit www.philips.com/respironics

Respironics and the names of Respironics products included here are trademarks of Respironics, Inc. and its affiliates.
Names of non-Respironics products are trademarks of their respective manufacturers.



© 2010 Koninklijke Philips Electronics N.V. All rights are reserved.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.
CAUTION: US federal law restricts these devices to sale by or on the order of a physician.

Hoech JJ 04/21/10 MCI 4103426 PN 1035985