

DeVilbiss Serenity[®] Nasal CPAP Mask

Simple and Intuitive Design

The simple and intuitive design of this economy nasal CPAP mask makes it the perfect choice. The Serenity mask features easy to size headgear and **interchangeable cushions and forehead pads**. The cushions are available in various materials and sizes. **The patented adjustable exhalation port** increases comfort by allowing the patient to select the direction the air exits the mask – no more blowing on the patient's chest or bed partner. **The lightweight mask is comfortable and provides a good seal.**



DeVilbiss Serenity
CPAP Mask
9352D

Performance Features

- Patented adjustable forehead rest takes pressure off the bridge of the nose for maximum patient comfort
- Patented adjustable vent allows patient to direct escaping air away from their chest or bed partner
- Snap together components are easy to disassemble and reassemble and minimizes likelihood of breakage
- Solid and reliable thermoplastic frame will hold up to repeated cleanings



Design Features

- 100% latex free
- Soft replaceable cushion provides a comfortable low pressure seal. Available in standard and shallow sizes and silicone or gel material
- Soft replaceable silicone or gel forehead pads
- Tubing runs over the head for improved tubing management and patient comfort
- Minimal headgear offers a secure fit and features a quick release buckle
- Supplemental oxygen port

Reimbursement information

| Product Description | Code Description | HCPCS Code | Frequency | Allowable Range |
|----------------------|--|------------|----------------|-----------------------------------|
| Serenity Mask | Nasal Application Device (mask or cannula type), used with positive airway pressure device, with or without head strap | A7034 | 1 per 3 months | \$99.00 - \$117.64 initial set-up |
| Replacement Cushion | Replacement cushion for Nasal Application Device, each | A7032 | 2 per 1 month | \$34.45 - \$40.53 |
| Replacement Headgear | Headgear | A7035 | 1 per 6 months | \$33.79 - \$39.75 |
| Chin Strap | Chinstrap | A7036 | 1 per 6 months | \$15.47 - \$18.20 |

Accessory and Replacement Parts

| DEVILBISS PART # | PRODUCT DESCRIPTION |
|------------------|---|
| 9352D | Serenity Mask with Standard Size Silicone Cushion |
| 9352G | Serenity Mask with Standard Size Gel Cushion |
| 9352S | Serenity Mask with Shallow Size Silicone Cushion |
| 9352GS | Serenity Mask with Shallow Size Gel Cushion |
| 9352D-601 | Replacement Standard Size Silicone Cushion |
| 9352G-601 | Replacement Standard Size Gel Cushion |
| 9352S-601 | Replacement Shallow Size Silicone Cushion |
| 9352GS-601 | Replacement Shallow Size Gel Cushion |
| 9352D-604 | Silicone Forehead Pads |
| 9352G-604 | Gel Forehead Pads |
| 9352D-619 | Headgear |



DeVilbiss Healthcare | 100 DeVilbiss Drive | Somerset, PA 15501 USA
 800.338.1988 | 814.443.4881 | www.DeVilbissHealthcare.com