

## SomnoPlus<sup>®</sup> Nasal Mask

### Freedom and Flexibility to Fit Your Lifestyle



SomnoPlus  
Nasal Mask  
(WM25820)



The SomnoPlus Nasal Mask offers ultimate comfort by combining **flexibility with exceptional fit**. The lightweight design allows for **freedom of movement** for even the most active sleeper. The **extremely quiet exhalation system** allows the patient and partner to sleep soundly and comfortably. Available in four sizes, SomnoPlus features integrated pressure measurement ports to allow convenient pressure measurement at the mask.

**Designed to adjust easily and precisely**, the SomnoPlus Nasal Mask provides a custom fit with two-way adjustable forehead support and breathable headgear. SomnoPlus is made to fit every patient therapy and lifestyle.

### Performance Features

- **Ball-and-socket joint** allows even the most active sleeper to move freely and comfortably
- Exceptionally quiet exhalation system allows the patient and bed partner to sleep soundly with **minimal noise**
- **Integrated pressure measurement ports** provide convenient pressure measurement at the mask
- **Small, lightweight design** allows for comfortable extended use and increases patient compliance

## Reimbursement Information

Product Description	Code Description	HCPCS Coding	Frequency	Allowable Range
Mask Frame	Nasal interface (mask or cannula type) used with positive airway pressure device, with or without head strap	A7034	1 per 3 months	\$99.00 – \$117.60, initial setup
Headgear	Headgear used with PAP device, each	A7035	1 per 6 months	\$30.57 – \$35.97
Cushion	Cushion for use on nasal mask interface, replacement only, each	A7032	2 per 1 month	\$34.45 – \$40.53

## Design Features

- Velcro headgear straps adjust easily to provide comfortable fit, seal and removal
- **Two-way adjustable forehead support** provides customized fit
- Broad forehead support minimizes uncomfortable pressure points on the forehead and nose by allowing even distribution of pressure
- Textured mask surface prevents sweating and moisture buildup
- Easy-to-clean mask is **dishwasher safe**; headgear is easily cleaned in the washing machine
- **Minimal parts** make mask easy to assemble
- Available in **four standard sizes** – XS, S, M and L
- Latex free

## Parts and Accessories

PRODUCT DESCRIPTION	DEVILBISS PART #
<b>Complete System</b>	
<b>Frame, Cushions and Headgear</b>	
SomnoPlus Nasal Mask, XS	WM25800
SomnoPlus Nasal Mask, S	WM25810
SomnoPlus Nasal Mask, M	WM25820
SomnoPlus Nasal Mask, L	WM25830
<b>Replacement Frame</b>	
<b>Frame Only (Mask without cushions or headgear)</b>	
SomnoPlus Frame Set, S	WM15709
SomnoPlus Frame Set, M	WM15710
SomnoPlus Frame Set, L	WM15711
<b>Replacement Cushions</b>	
EasyFit / SomnoPlus Nasal Cushion, Sil, S	DV15812
EasyFit / SomnoPlus Nasal Cushion, Sil, M	DV15813
EasyFit / SomnoPlus Nasal Cushion, Sil, L	DV15814
<b>Replacement Headgear</b>	
SomnoPlus Headgear	WM25150



DeVilbiss Healthcare LLC | 100 DeVilbiss Drive | Somerset, PA 15501 USA  
 800.338.1988 | 814.443.4881 | [www.DeVilbissHealthcare.com](http://www.DeVilbissHealthcare.com)