



**PHILIPS**  
RESPIRONICS

*DreamStation*

CPAP

## DreamStation helps support long-term patient use while creating **new efficiencies for your business.**

**Connectivity** options offer more monitoring choices than any leading CPAP brand with integrated Bluetooth® on every device, and offering optional Wi-Fi and cellular modems.\*

**Performance Check** simplifies in-home device evaluation and can help reduce the number of normally functioning devices that are returned for service.

**Daily Progress Feedback** displays a simple trend of patients' nightly use to help track progress on the path to adherence.

### **DreamStation features**

- Advanced event detection
- Mask type resistance control
- Mask fit check
- Oximetry capable
- Encore patient management system with optional Adherence Profiler
- DreamMapper patient self-management system compatibility

### **Humidification**

- Fixed, adaptive, heated tube
- Pre-heat humidifier option
- Modular design with 1 piece, dishwasher-safe water chamber

### **Protection programs**

- Flex Promise program
- Bi-level Rescue program
- "One-Touch" warranty programs

### **CPAP features**

- C-Flex
- SmartRamp
- Basic AHI/Leak reporting
- Standard integrated Bluetooth

Provided for your reference by

**cpapXchange.com**

\*Internal assessment of 2015 competitive CPAP data comparing to ResMed Aircurve10/Aircurve10 platform and Fisher & Paykel Icon series platform. The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and the use of such marks by Philips Respironics is under license. Other trademarks and trade names are those of their respective owners.

# DreamStation CPAP

Comfort features			Encore data capture and reporting	
Modes	Features	Flex Comfort	Advanced Event Detection	Flow Waveforms
Fixed CPAP	N/A	C-Flex	N/A	N/A

## Specifications

### DreamStation CPAP

Pressure Range	4 to 20 cm H2O
Flex pressure relief	0 to 3
Ramp time	0 to 45 min (5-minute increments)
Starting ramp pressure	4 cm H2O to CPAP Pressure (CPAP)
Humidification	Heated humidification: Fixed, Adaptive, Heated Tube modes; Heated tube (optional)
Data storage capacity (minimum)	SD card: > 1 year On-board: <6 months
Altitude compensation	Automatic
Filters	Reusable pollen; Disposable Ultra-fine
Device controls	LCD, control dial/push button
Warranty	2 years (US)
Dimensions	15.7x19.3x8.4cm (w/o humidifier) 29.7x19.3x8.4cm (w/ humidifier)
Weight	1.33kg/2.94lbs (w/o humidifier) 1.98kg/4.37lbs (w/humidifier) Incl. power supply
Electrical requirements	100 – 240 VAC, 50/60 Hz, 2.0-1.0 A

## Ordering information

DreamStation CPAP	Part number
DreamStation CPAP	DSX200S11
DreamStation CPAP with humidifier	DSX200H11
DreamStation CPAP with humidifier and heated tube	DSX200T11
Power supply	Part number
DreamStation 80-watt power supply	1118499
Power cord	1005894

## Ordering information

Humidification	Part number
DreamStation humidifier	DSXH
DreamStation water chamber	1122520
RP-DS Humidifier Dry Box Assembly	1120668
RP-DS Humidifier Flip Lid Seal	1120617
RP-DS Humidifier Dry Box Inlet Seal	1120613
Tubes	Part number
15MM Std Tube – DreamStation - RP	PR15
15MM Heated Tube	HT15
Filters	Part number
Pollen filter, reusable (1 per pack)	1122446
Ultra-fine filter, disposable (1 per pack)	1122447
Ultra-fine filter, disposable (2 per pack)	1122518
Ultra-fine filter, disposable (6 per pack)	1122519
Accessories	Part number
Link module	1120293
Nonin SpO2, assembly	1121694
Reuseable finger sensor (clip)	936
Oximetry finger sensor, adult, flex	953
Flexiwrap, sens tape, adult, 25/pk	954A
Reuseable finger sensor (soft), medium	1062050
Reuseable finger sensor (soft), small	1062051
Modems	Part number
Cellular modem, oximetry capable	100600C
Cellular modem, Non-oximetry	100605C
Wi-Fi accessory, oximetry capable	100700W
Wi-Fi accessory, Non-oximetry	100705W

© 2015 Koninklijke Philips N.V. All rights reserved. Specifications are subject to change without notice. Trademarks are the property of Koninklijke Philips N.V. (Royal Philips) or their respective owners.

www.philips.com



Caution: US federal law restricts these devices to sale by or on the order of a physician.

edoc EB 10/05/15 MCI 4106903  
1010 Murry Ridge Lane, Murrysville, PA 15668  
800-345-6443 • 724-387-4000