

## DreamStation Sleep Therapy System Configurations:

|  | DreamStation CPAP  | DreamStation CPAP Pro                                    | DreamStation Auto CPAP                                   | DreamStation BiPAP Pro                                   | DreamStation Auto BiPAP                                  |
|--|--|--|--|--|--|
| <b>Modes</b>                           | Fixed CPAP   | Fixed CPAP<br>CPAP-Check                                 | Fixed CPAP<br>CPAP-Check<br>Auto-CPAP                    | Fixed CPAP<br>Fixed BiPAP                                | Fixed CPAP<br>Fixed BiPAP<br>Auto-BiPAP                  |
| <b>Features</b>                        | n/a  | EZ-Start<br>Auto-Trial                                   | EZ-Start<br>Auto-Trial<br>Opti-Start                     | EZ-Start   | EZ-Start   |
| <b>Pressure Range</b>                  | 4 to 20 cm H <sub>2</sub> O                              | 4 to 20 cm H <sub>2</sub> O                              | 4 to 20 cm H <sub>2</sub> O                              | 4 to 25 cm H <sub>2</sub> O                              | 4 to 25 cm H <sub>2</sub> O                              |
| <b>Flex Comfort</b>                    | C-Flex   | C-Flex<br>C-Flex+  | C-Flex<br>C-Flex+<br>A-Flex                              | C-Flex<br>C-Flex+<br>Bi-Flex                             | C-Flex<br>C-Flex+<br>Bi-Flex                             |
| <b>Humidification Modes (3)</b>        | Fixed<br>Adaptive<br>Heated Tube                         | Fixed<br>Adaptive<br>Heated Tube                         | Fixed<br>Adaptive<br>Heated Tube                         | Fixed<br>Adaptive<br>Heated Tube                         | Fixed<br>Adaptive<br>Heated Tube                         |
| <b>Mask Type Resistance Control</b>    | Standard   | Standard   | Standard   | Standard   | Standard   |
| <b>Ramp (0-45 minutes)</b>             | Ramp<br>SmartRamp  | Ramp<br>SmartRamp  | Ramp<br>SmartRamp  | Ramp<br>SmartRamp  | Ramp<br>SmartRamp  |
| <b>Altitude Compensation</b>           | Automatic  | Automatic  | Automatic  | Automatic  | Automatic  |
| <b>Advanced Event Detection</b>        | n/a  | Standard   | Standard   | Standard   | Standard   |
| <b>Flow Waveforms</b>                  | n/a  | Standard   | Standard   | Standard   | Standard   |
| <b>Modem Compatibility</b>             | Standard<br>(optional Cellular Modem or Wi-Fi Accessory) | Standard<br>(optional Cellular Modem or Wi-Fi Accessory) | Standard<br>(optional Cellular Modem or Wi-Fi Accessory) | Standard<br>(optional Cellular Modem or Wi-Fi Accessory) | Standard<br>(optional Cellular Modem or Wi-Fi Accessory) |
| <b>Built-in Bluetooth Connectivity</b> | Standard   | Standard   | Standard   | Standard   | Standard   |
| <b>DreamMapper Compatibility</b>       | Standard<br>(by mobile app or web)                       | Standard<br>(by mobile app or web)                       | Standard<br>(by mobile app or web)                       | Standard<br>(by mobile app or web)                       | Standard<br>(by mobile app or web)                       |
| <b>Oximetry Compatibility</b>          | Standard<br>(optional Nonin accessories)                 | Standard<br>(optional Nonin accessories)                 | Standard<br>(optional Nonin accessories)                 | Standard<br>(optional Nonin accessories)                 | Standard<br>(optional Nonin accessories)                 |
| <b>Check Mask Fit feature</b>          | Standard   | Standard   | Standard   | Standard   | Standard   |
| <b>Performance Check</b>               | Standard   | Standard   | Standard   | Standard   | Standard   |