

DreamStation Sleep Therapy System Configurations:

	DreamStation CPAP	DreamStation CPAP Pro	DreamStation Auto CPAP	DreamStation BiPAP Pro	DreamStation Auto BiPAP
Modes	Fixed CPAP	Fixed CPAP CPAP-Check	Fixed CPAP CPAP-Check Auto-CPAP	Fixed CPAP Fixed BiPAP	Fixed CPAP Fixed BiPAP Auto-BiPAP
Features	n/a	EZ-Start Auto-Trial	EZ-Start Auto-Trial Opti-Start	EZ-Start	EZ-Start
Pressure Range	4 to 20 cm H ₂ O	4 to 20 cm H ₂ O	4 to 20 cm H ₂ O	4 to 25 cm H ₂ O	4 to 25 cm H ₂ O
Flex Comfort	C-Flex	C-Flex C-Flex+	C-Flex C-Flex+ A-Flex	C-Flex C-Flex+ Bi-Flex	C-Flex C-Flex+ Bi-Flex
Humidification Modes (3)	Fixed Adaptive Heated Tube	Fixed Adaptive Heated Tube	Fixed Adaptive Heated Tube	Fixed Adaptive Heated Tube	Fixed Adaptive Heated Tube
Mask Type Resistance Control	Standard	Standard	Standard	Standard	Standard
Ramp (0-45 minutes)	Ramp SmartRamp	Ramp SmartRamp	Ramp SmartRamp	Ramp SmartRamp	Ramp SmartRamp
Altitude Compensation	Automatic	Automatic	Automatic	Automatic	Automatic
Advanced Event Detection	n/a	Standard	Standard	Standard	Standard
Flow Waveforms	n/a	Standard	Standard	Standard	Standard
Modem Compatibility	Standard (optional Cellular Modem or Wi-Fi Accessory)	Standard (optional Cellular Modem or Wi-Fi Accessory)	Standard (optional Cellular Modem or Wi-Fi Accessory)	Standard (optional Cellular Modem or Wi-Fi Accessory)	Standard (optional Cellular Modem or Wi-Fi Accessory)
Built-in Bluetooth Connectivity	Standard	Standard	Standard	Standard	Standard
DreamMapper Compatibility	Standard (by mobile app or web)	Standard (by mobile app or web)	Standard (by mobile app or web)	Standard (by mobile app or web)	Standard (by mobile app or web)
Oximetry Compatibility	Standard (optional Nonin accessories)	Standard (optional Nonin accessories)	Standard (optional Nonin accessories)	Standard (optional Nonin accessories)	Standard (optional Nonin accessories)
Check Mask Fit feature	Standard	Standard	Standard	Standard	Standard
Performance Check	Standard	Standard	Standard	Standard	Standard