

Provided for your reference by  
**cpapXchange.com**

# Economies of scale

**PROBASICS**™  
by PMI  
**Zzz-PAP**™



# Introducing the new **Zzz-PAP™** CONTINUOUS POSITIVE AIRWAY PRESSURE UNIT

The Probasics Zzz-PAP™ by PMI delivers performance, simplicity, style, and value in a very small package. The Zzz-PAP offers easy patient setup with a 4-18 cmH<sub>2</sub>O pressure range, a comfortable 0-45 minute ramp feature and manual altitude adjustment. PMI also offers the optional Probasics Zzz-PAP Heated Humidifier that quickly and easily integrates to the Zzz-PAP unit providing excellent relief from the symptoms of dryness and congestion associated with CPAP therapy. Weighing less than two pounds and coming standard with an auto-adjusting 100 Volt/240 Volt Universal Power Supply, the compact Zzz-PAP is not only an excellent choice for your bedside but it's perfect for travel as well. Best of all, the Probasics Zzz-PAP is backed by a two-year manufacturer's warranty which we all know is the most necessary feature for providing both the patient and the dealer with peace of mind.



7501

## **Zzz-PAP™** INTEGRATED HEATED HUMIDIFIER

Designed for use with the Probasics Zzz-PAP™ Mini-CPAP unit, PMI offers an Integrated Heated Humidifier system to relieve nasal irritation and dryness which result from constant airflow by supplying the required humidity to the upper airway. With a reliable and safe design which features overheat protection, the Probasics Zzz-PAP Heated Humidifier can drastically improve the patient's comfort level and improve patient compliance.



7600



7500

---

7500	Zzz-PAP CPAP Unit with Interated Heated Humidifier
7501	Zzz-PAP CPAP Unit alone
7601	Zzz-PAP Interated Heated Humidifier alone

---

# SPECIFICATIONS



## Continuous Positive Airway Pressure

### Easy-to-carry design

It is compact and easy to carry with its small size and light weight...Perfect for travel.

### LCD display with back light and easy-to-use interface

The LCD display is easy to read, and the intuitive user interface is very simple to use.

### Easy Pressure level adjustment

Operating pressure, ranging from 4 to 18 cm H<sub>2</sub>O, can be easily and securely set to comply with the physician's prescription while in provider mode.

### Easy Ramp time setting and adjustment

Deluxe ramp capability allows a 0 - 45 minute adjustable time period to help user fall to sleep easier.

### Compliance meter

The number of patient use hours can be easily accessed for reimbursement or therapeutic recording purposes.

### Extremely quiet operation

Performs at just 30 dBA at 10.0 cmH<sub>2</sub>O.

Pressure Range	4.0 - 18.0 cmH <sub>2</sub> O (0.5 cmH <sub>2</sub> O increment)
Pressure Accuracy	+2 -1 cmH <sub>2</sub> O
Ramp Time	0 - 45 min (5 minute increments)
Altitude Compensation	Level 1 ~ 8 for 0 ~ 8000 ft. (manual setting)
Device Operating Meter	9999.9 hr
Heated Humidifier	Yes
Sound Level	30 dBA @ 10 cmH <sub>2</sub> O, one meter distance
Dimension (W x D x H)	14.5 x 13.0 x 10.0 cm 5.7" x 5.1" x 3.9"
Weight	800 g (estimation)
Power Supply	100V-240V, 50/60Hz, Universal Power Supply
Operating Temperature	5°C to 35°C / 41°F to 95°F
Environment Humidity	15% to 95% RH



## Integrated Heated Humidifier

- Integrates with Probasics Zzz-PAP mini CPAP unit
- Easy to use rotary knob for heating level adjustment
- Dishwasher safe water chamber has an easy to remove base plate for cleaning
- Safe & reliable with overheat protection
- Indicators clearly show when Power is ON and Heater is ON

Power Supply	120VAC, 50/60Hz    230VAC, 50/60Hz
Water Capacity	450 ml
Heater Settings	1 - 6 (Heater plate temperature 86°F to 140°F / 30°C to 60°C)
Warm-Up Time	< 30 minutes
Pressure Drop	0.2 cmH <sub>2</sub> O @ 60LPM
Humidity Output	10 - 40 mgH <sub>2</sub> O/L
Indicator	Power On indicator, Heater On indicator
Dimension (W x D x H)	Heater: 6.8" x 10" x 1.6" (17.2 x 25.5 x 4 cm) Water Chamber: 4.5" x 5.9" x 4.2" (11.5 x 14.9 x 10.6 cm)
Weight	Heater: 380 g (0.84 lb) Water Chamber: 240 g (0.53 lb)

All specifications are subject to change without notice.

0607



Probasics™

15 South Main Street

Second Floor, Unit F

Marlboro, NJ 07746

Tel. 800.438.0503

Fax 732.683.1945

www.pmiusa.biz