

**ADAM™**

## INTERFACE SYSTEM

**Proven Interface Technology.**

Lighter and more comfortable than traditional nasal mask interfaces, the *ADAM™* Interface System (Airway Delivery and Management) uses Nasal Pillows, offering patients a unique alternative for noninvasive ventilation therapy. The Nasal Pillows gently seal at the patient's nares to reduce leaks. The *ADAM* interface includes Velcro® fasteners that permit maximum versatility in fitting and adjustment.



*ADAM*  
Interface  
with  
*KnightWear®*  
Headgear



*ADAM*  
Interface  
with  
*CoolKnight™*  
Headgear



In addition to the *ADAM* system, Puritan Bennett offers a range of interface options, including the innovative *Breeze SleepGear™* Interface System shown here.

**Nasal Pillows minimize leaks and give CPAP patients a comfortable alternative to traditional masks.**

Nasal Pillows come in seven sizes to give your patients a wide choice in fit.

**Swivel adapters relieve stress between the tubing and shell to help keep pillows in place.**

Tubing routes over patient's head for greater freedom of movement.

**Four lightweight, adjustable headgear options are available for sleep lab and home use.**

A calibrated exhalation hole in the shell allows carbon dioxide venting.

# ADAM

## INTERFACE SYSTEM

### Ordering Information

A complete patient circuit consists of a circuit assembly, headgear and nasal pillows. Other accessories can be ordered as needed. Circuit assemblies include 6' smooth bore tube, swivel adapter and shell. Hose attachment fits 22-mm OD connection.

Description	Catalog Number
<i>ADAM</i> narrow/large <i>SnugFit</i> kit, includes narrow shell, large pillows, <i>SnugFit</i> Headgear, flex tube, swivel adapter, angle adapter and 22-mm connector	Y-100874-00
<i>ADAM</i> narrow/medium <i>KnightWear</i> kit, includes narrow shell, medium pillows, <i>KnightWear</i> Headgear, flex tube, swivel adapter, angle adapter and 22-mm connector	Y-100873-00
<i>ADAM</i> standard/medium <i>SnugFit</i> kit, includes standard shell, medium pillows, <i>SnugFit</i> Headgear, flex tube, swivel adapter, angle adapter and 22-mm connector	Y-100871-00
<i>ADAM</i> standard/large <i>SnugFit</i> kit, includes standard shell, large pillows, <i>SnugFit</i> Headgear, flex tube, swivel adapter, angle adapter and 22-mm connector	Y-100872-00
Nasal ventilation circuit, includes standard shell, swivel adapter, short connecting hose	S-133283-00
<i>CoolKnight</i> Headgear with Breath-O-Prene	Y-101564-00
<i>KnightWear</i> Cotton Headgear	Y-100028-00D
Reusable <i>SnugFit</i> Headgear, one size fits all, with shell strap, black	Y-133317-00E
Disposable <i>SnugFit</i> Headgear, one size fits all, with shell strap, white (5/pkg)	S-133321-00
Circuit assembly without pressure tap, for home use	S-133293-00
Circuit assembly with pressure tap, for sleep lab use	S-133294-00
Velcro patches for <i>SnugFit</i> Headgear (6/pkg)	S-290560-00A
Nasal Pillows, small, one pair, blue	S-616323-00A
Nasal Pillows, medium, one pair, clear	S-616324-00C
Nasal Pillows, large, one pair, green	S-616325-00A
Nasal Pillows, extra large, one pair, lavender	Y-101969-00
Dilator Nasal Pillows, small, one pair, grey	S-616355-00C
Dilator Nasal Pillows, medium, one pair, white	S-616356-00C
Dilator Nasal Pillows, large, one pair, yellow	S-616357-00C

### Accessories and Replacement Parts

CPAP shell with thermistor ports	S-231699-00B
CPAP calibration shell	S-231702-00
CPAP shell, narrow	S-232101-00
Narrow clinical shell with thermistor ports	S-232102-00
CPAP shell without pressure tap	S-231700-00B
CPAP shell with pressure tap	S-231701-00
Vent shell, standard	S-231703-00
Vent shell, narrow	S-232100-00
Pillow sizing gauge (10/pkg)	S-616468-PG
Oxygen adapter	S-616155-03B
Swivel adapter	S-616329-00
Angle adapter	S-616369-00
Chin strap, black	S-133318-00
Shell strap, black	S-133319-00
Smooth bore air supply tubing, 6'	S-133292-00
Smooth bore air supply tubing, 8'	S-776274-00C
Air supply tubing, 15", with swivel adapter	S-133324-00
Flex tube assembly with swivel adapter	S-133350-00
Disposable smooth bore air supply tubing, 6', 6/cs	S-133325-00
<i>ADAM</i> CPAP instruction sheet	Y-799702-00E

**tyco**

Healthcare

**Puritan  
Bennett**

4280 Hacienda Drive  
Pleasanton, CA 94588  
Tel 925.463.4000  
Toll Free 1.800.635.5267  
www.puritanbennett.com

Tyco Healthcare UK LTD.  
154 Fareham Road  
Gosport, UK PO13 0AS  
Tel +44.1329.224000

Caution: Federal (U.S.A.)  
law restricts this device  
to sale by or on the order  
of a physician.

Specifications are subject  
to change without notice.

© 2003 Nellcor Puritan  
Bennett Inc. All rights  
reserved.

A.z 1039-0203  
ST02400