



## Classic comfort

**ComfortClassic nasal mask** provides quality performance in a no-nonsense package.


Sometimes a simple solution is the best. No wonder this uncomplicated mask remains a favorite to those looking for a good seal, reduced complexity, and great value.

ComfortClassic sets a standard of comfort and convenience for virtually any of your male or female patients. With a dual-layer silicone cushion, the outer membrane creates a comfortable seal, while the inner layer provides stability. A dual-position gel spacer allows the mask to adjust to the patient's profile. It also conforms to the wearer's forehead for even greater comfort and stability. Round out the package with our deluxe headgear, and you have a high-performance, affordable mask that helps you fit 100% of your patients, 100% of the time... while also helping your bottom line.

Provided for your reference by  
**cpapXchange.com**

**PHILIPS**

**RESPIRONICS**



ComfortClassic features our one-size-fits-most deluxe headgear for increased stability and improved seal. With four adjustment points, it distributes tension evenly for optimum comfort. The headgear can be used with most of our masks.



Item				
<b>ComfortClassic mask</b>	<b>Small</b>	<b>HCPCS code</b>	<b>Medium</b>	<b>HCPCS code</b>
With headgear	1007963	A7034,A7035	1007964	A7034,A7035
Mask only	1007966	A7034	1007965	A7034
<b>Parts and accessories</b>			<b>Part number</b>	<b>HCPCS code</b>
Deluxe headgear			1002800	A7035
Deluxe headgear, small			1018755	A7035
Gel forehead spacer kit (standard and short)			1002464	N/A
Split washer kit (5 per kit)			1002720	N/A
Exhalation elbow kit			1004395	A7045
Foam forehead spacer (5 per pack)			1008651	N/A

Respironics and ComfortClassic are trademarks of Respironics, Inc. and its affiliates.



©2009 Koninklijke Philips Electronics N.V. All rights are reserved. Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication. CAUTION: US federal law restricts these devices to sale by or on the order of a physician.

Philips Healthcare is part of Royal Philips Electronics  
 Europe, Middle East, Africa:  
 +33 1 47 52 30 00  
 Asia Pacific:  
 +852 3194 2280

Philips Respironics  
 1010 Murry Ridge Lane  
 Murrysville, PA 15668  
 +1 724 387 4000  
 +1 800 345 6443