



**Note:** For proper mask sizing print out this page at 100%.

Do NOT select the "Shrink oversized pages to paper size" option.

The box in the upper right corner should measure 1 inch square.



The diagram shows two teardrop-shaped mask templates labeled 'M' (Medium) and 'S' (Small). Each template has a horizontal line across its widest part with arrows pointing outwards and the text 'cut out along lines' below it. To the left of the 'M' template are two measurement diagrams: the top one shows a face with 'L' for vertical distance from the bridge of the nose to the upper lip and 'W' for the width of the upper lip; the bottom one shows a mask with 'L' for vertical distance from the top edge to the bottom edge and 'W' for the width at the bottom edge.

ComfortFusion nasal mask sizing gauge  
www.philips.com/comfortfusion

PHILIPS  
RESPIRONICS

Provided for your reference by  
**cpapXchange.com**