

1 inch



● Full Face Delivery ● Direct Nasal Delivery  
● Nasal Delivery ● Oral Delivery

L  
M  
S  
XS

**Infinity<sup>481</sup>**  
DIRECT NASAL MASK



**FlexiFit<sup>407</sup>**  
PREMIUM NASAL MASK



**FlexiFit<sup>406</sup>**  
PETITE NASAL MASK

L  
M  
S

Provided for your reference by  
**cpapXchange.com**

### Fisher & Paykel Multi-Mask Sizing Gauge

1. Print document - Note do not choose to scale document to paper size when printing. (Square at top of page should be one inch wide/tall)
2. Cut out paper inside of mask outlines
3. Choose correct size according to fit (406/407) or size lines (431/432/481).

**FlexiFit<sup>431/432</sup>**  
FULL FACE MASK

**SleepStyle<sup>™</sup>**  
SLEEP PERFORMANCE FOR  
AN ENERGIZED LIFESTYLE

**Fisher & Paykel**  
HEALTHCARE