

Technology Highlight: **Frame**

“Extremely light and you don’t feel claustrophobic.”


“Softer on the forehead; no marks on the nose bridge.”

Reference: ResMed user trials



Mirage™ FX is our most streamlined, unobtrusive nasal mask on the market; users will love the feeling of openness and clear line of sight.

The lightweight, durable frame is strong enough to stabilise the mask on the face, yet soft enough to accommodate user movement during sleep.



Flexible, winged forehead support

allows a dynamic range of fit and adjustment; the plastic does not come into contact with the face.

Mask weight: 70g
(Total mass including headgear)



Over-moulded headgear loops

simplify attachment and detachment of headgear, and retain optimal headgear settings.

The loops are integrated into the streamlined frame design, eliminating the need for extra parts, that is, headgear clips.



Effortless and
effective

Mirage™ FX
NASAL MASK

www.resmedmiragefx.com