





Ultra Mirage™ II
NASAL MASK

Mirage Activa™
NASAL MASK


Provided for your reference by
cpapXchange.com




↑
STANDARD




↑
LARGE



↑
SHALLOW



↑
SHALLOW WIDE
(Ultra Mirage II only)



Actual Size: 5 1/2" (Width) x 7" (Height)

Ruler is for scaling purposes.

How to use the fitting template:

1. Cut out the template along the dotted line
2. Cut out nasal holes
(hint: fold each nasal hole in half to cut)
3. Follow directions on template to determine appropriate cushion size

RESMED