

# Nasal Aire II<sup>®</sup> Headgear

## Multiple Wearing Options

### Standard Method :

All components included in Package



Use headband , two (2) tubing holder rings , neckstrap and bolo .

### Florida Method :

\* Optional headgear



Use Headband , two (2) figure 8 loops and bolo.

### Ponytail Method :

All components included in Package



Use Headband , two (2) tubing holder rings and bolo.

### Breatheband Method :

\* Optional headgear



Use Breatheband and neckstrap.

**Chinstrap Method :**

\* Optional headgear



Use chinstrap and neckstrap

**Figure 8 Method :**

All components packaged



**Use figure 8 loop**

\* Denotes optional item. Sold Separately.