



# Customized to Energize.

## Zest™ Range

### NASAL MASK

Compact, quiet, and effortless to fit and use, the Zest™ Nasal Mask Range offers superb comfort. The innovative design of the Zest™ Nasal Masks combine the Easy-Clip Silicone Seal and Advanced Air Diffuser with proven technology that makes Fisher & Paykel Healthcare masks perform so well.

Simple to both use and seal, the Zest™ Range of Nasal Masks is your complete solution for patient satisfaction!

***With Zest, the rest is easy***

Easy-Clip™ Silicone Seal  
with VTS™



**Zest™** Petite




**Zest™**




**Zest™** Plus

# Nasal Delivery – Zest™



The Zest™ Range of Nasal Masks offers an over-the-nose mask solution that fits better, feels lighter and is easier to use. The mask includes more features that make Fisher & Paykel Healthcare masks so popular. The Zest™ Nasal Mask comes in three separate sizes, as an Out-of-Box Solution fully assembled and ready for use.

## UNIQUE MASK FEATURES FOR GREATER FREEDOM, COMFORT AND COMPLIANCE




### FlexiFit™ Technology Auto-contouring fit

With Foam Cushion FlexiFit™ Technology, automatically adjusts to a wide range of nasal and facial contours with simple one-step ease of fitting and no complicated adjustments.




### FlexiFoam™ Cushion Active-contouring cushion

The FlexiFoam™ Cushion provides greater comfort with the unique design of the soft, light and supportive foam. The FlexiFoam™ Cushion features the advanced active contouring that conforms naturally to your patient's face.




### Glider™ Strap Anti-leak system

The sliding Glider™ Strap provides freedom of movement while maintaining a seal. Mask force remains equal when the head rotates from side to side, minimizing the occurrence of leaks.



### Easy-Clip Silicone Seal with VTS™ Effortless performance

The Easy-Clip Silicone Seal allows one-step connection or removal for quick and effortless cleaning. The Variable Thickness Silicone (VTS™) provides an ultra-soft fit for improved sealing and comfort.




### Advanced Air Diffuser Quiet usage

The maintenance-free Air Diffuser disperses the air-flow away from the sleeper's partner, and minimizes exhaust flow noise.




### Stretchgear™ Headgear Deluxe mask suspension system

A comfortable and breathable elastic headgear, specially designed to maintain a firm seal.



### Glider Clip / Hidden Strap Release One-step usage

Easy one-step removal of headgear with Glider Clip or Hidden Strap Release. Enables reassembly without having to make fitting adjustments.



\*According to the Fisher & Paykel Healthcare Anthropometric Database.

## CONTENTS

- Complete mask with headgear assembled and ready for use
- Three sizes: Petite or Standard or Plus
- User-friendly Fitting Guide

**REORDER CODE:**  
400439 Zest™ Petite  
400440 Zest™  
400441 Zest™ Plus

FOR MORE INFORMATION, PLEASE CONTACT YOUR LOCAL FISHER & PAYKEL HEALTHCARE REPRESENTATIVE:

**Manufacturer** Fisher & Paykel Healthcare Ltd, 15 Maurice Paykel Place, East Tamaki, Auckland 2013, PO Box 14 348, Panmure, Auckland 1741, New Zealand  
Tel: +64 9 574 0100 Fax: +64 9 574 0158 Email: info@fphcare.com Web: www.fphcare.com

**Australia** Tel: +61 3 9879 5022 Fax: +61 3 9879 5232 **Austria** Tel: 0800 29 31 23 Fax: 0800 29 31 22 **Benelux** Tel: +31 40 216 3555 Fax: +31 40 216 3554 **China** Tel: +86 20 3205 3486 Fax: +86 20 3205 2132  
**France** Tel: +33 1 6446 5201 Fax: +33 1 6446 5221 **Germany** Tel: +49 7182 93777 0 Fax: +49 7182 93777 99 **India** Tel: +91 80 4123 6041 Fax: +91 80 4123 6044 **Irish Republic** Tel: 1800 409 011  
**Italy** Tel: +39 06 7839 2939 Fax: +39 06 7814 7709 **Japan** Tel: +81 3 3661 7205 Fax: +81 3 3661 7206 **Middle East** Tel: +9 71 4 360 0849 Fax: +9 71 4 368 8110 **Northern Ireland** Tel: 0800 132 189  
**Spain** Tel: +34 902 013 346 Fax: +34 902 013 379 **Sweden** Tel: +46 8 564 76 680 Fax: +46 8 36 63 10 **Switzerland** Tel: 0800 83 47 63 Fax: 0800 83 47 54 **UK (EU Authorised Representative)**  
Fisher & Paykel Healthcare UK Ltd, Unit 16, Cordwallis Park, Clivemont Road, Maidenhead, Berkshire SL6 7BU, UK Tel: +44 1628 626 136 Fax: +44 1628 626 146 **USA** Tel: +1 800 446 3908 or +1 949 453 4000  
Fax: +1 949 453 4001

© Fisher & Paykel Healthcare Limited